ozū

MANUAL DE INSTRUCCIONES

# ELECTROESTIMULADOR 3 EN 1 TENS EMS FITNESS

UM-004

WELL-BEING SERIES

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# ozū

USER MANUAL

## 3 IN 1 TENS EMS FITNESS ELECTROSTIMULATOR

UM-004

WELL-BEING SERIES

## IMPORTANT SAFETY INSTRUCTIONS

Dear customer, we thank you deeply for your purchase of this device. Please read this User Manual carefully and keep it for further reference. When using the product, all instructions and basic precautions should always be followed, including the following:

### DANGER

To reduce the risk of electric shock:

Always unplug this product from the electrical outlet immediately after use and before cleaning.

### WARNINGS:

To reduce the risk of burns, fire, electric shock, or injuries:

1. The product should never be left unattended when plugged in. Unplug from outlet when the device is not use and before switching in or out parts.
2. Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
3. Do not carry this product by its supply cord or use the cord as a handle.
4. Never operate the product with the air openings blocked. Keep the air openings free of lint, hair, or similar.
5. Do not drop or insert any object into any opening.
6. Do not use outdoors.
7. Do not operate where aerosol (spray) products are being used or where oxygen is administered.
8. Be sure to pull the plug out of the socket when you reach in. And don't pull it out when your hands are wet. That could cause electric shock and injury.
9. Remove the power plug from the socket immediately after a power failure. Otherwise, it may cause electric shock, short circuit and fire.
10. Keep the cord away from heated surfaces. If a fault is found, stop using immediately and cut off the power supply. Avoid accidents or injuries. Please do not repair it by yourself, please contact the seller or manufacturer in time.
11. The appliance must only be supplied at safe low voltage corresponding to the marking on the appliance.
12. The appliance is only to be used with the power supply unit provided with the appliance.
13. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
14. Special attention must be paid to components most susceptible to wear and tear. Replace defective components immediately and/or keep the equipment out of use until repair.
15. Use this appliance only for its intended use as described in this manual.
16. Children should be supervised to ensure that they do not play with the appliance.
17. Cleaning and user maintenance shall not be made by children without supervision.

\*If you feel any abnormality during use, please stop using immediately and consult

your doctor. Otherwise may it be harmful to health.

The following must not use this product:

1. People who are susceptible to electromagnetic malfunctions such as cardiac pacemakers and implantable medical electronic devices.
2. People who are susceptible to lose their sense of temperature.
3. People with malignant tumors.
4. People with heart disease.
5. A person who is early in pregnancy or shortly after giving birth.
6. A person who has an abnormal or curved spine.
7. A person with a disturbance of perception caused by a high degree of failure in the peripheral circulation in diabetes.
8. Someone with a scar on their skin.
9. People who have been medically prescribed to rest and people who are in poor health.
10. People with osteoporosis, severe fractures, sprain, and muscle pain, etc., acute headache.
11. People with 38 °C temperature (fever).
12. Acute inflammatory symptoms, (tiredness, chills, changes in blood pressure, etc.) or going through a time of intense, debilitating stress.
13. People with tendonitis.
14. Others treated at medical facilities.

## Precautions

1. This product is a household electronic personal care product, not a substitute for medical equipment, and not suitable for diagnosis, treatment and other commercial purposes.
2. Before use, please remove the accessories on the massage part or other items that may affect the massage.
3. Please stop using if there is pain or discomfort in the massage area when using the product.
4. When using the product, if there is any abnormality in the product, please turn off the product immediately, cut off the power supply and contact customer service, and do not disassemble and repair it by yourself.
5. When using in other specific environments, please comply with safety regulations.
6. Please use the adapter that matches the input voltage of the product to connect the product, and pay attention to the safety of the input electricity. It is recommended not to connect at the same outlet simultaneously with other devices.
7. Please keep this product away from water, fire and corrosive environments.
8. Do not use this product in one place for a long time to prevent skin damage and muscle pain caused by excessive massage.
9. Do not use this product while focusing on precise operations or while driving.
10. Do not use it if the massage part has undergone surgery, or if there is discomfort such as skin damage, redness, inflammation, congestion, etc.
11. Patients with heart disease or medical history, acute disease, skin disease, malignant tumor, blood disease such as bleeding tendency, osteoporosis, soft

tissue disease, stroke and other physical abnormalities should not use this product.

12. People with slow movement, physical and mental discomfort or intellectual disabilities should not use it. If you need to use it, please follow your doctor's advice.
13. It is not recommended to use the massager during pregnancy and menstruation. If you need to use it, please follow your doctor's advice.
14. When storing the product, please place it out of the reach of children. Children under the age of 6 are prohibited to use, and children aged 6-12 must use it under the guidance and supervision of their parents or persons with the ability to supervise.

## 1. Product Introduction

- The Transcutaneous Electrical Nerve Stimulators is a portable and battery powered multifunctional device.
- The Transcutaneous Electrical Nerve Stimulators (model: UM-004) offering Transcutaneous Electronic Nerve Stimulator (TENS), Electronic Muscle Stimulator (EMS) and FITNESS. It has 50 operation modes, which can give certain electrical pulse through electrode pads placed on the skin to help users to enjoy body stimulation.
- The electronic stimulatory module has the operating elements of ON/OFF Key, Display screen, Mode Selection key and Intensity Modification keys.
- The LCD display screen can show selected mode and program, output intensity, stimulate frequency, and time remaining of an application mode.
- The device is equipped with accessories of electrode pads, electrode wires, and batteries. The electrode wire is used to connect the patches to the main unit.
- The electrode pads are complying with the biocompatibility standards ISO 10993-5 (Cytotoxicity) and ISO 10993-10 (Irritation and Sensitization), are interchangeable.

## 2. Product Intended use

- TENS (Mode 1~20)  
To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, arm, leg, and foot due to strain from exercise or normal household work activities by applying current to stimulate nerve.
- EMS (Mode 21~40)  
It is intended for relaxation of muscle spasm, prevention or retardation of disuse atrophy, increase local blood circulation, and muscle re-education.
- FITNESS (Mode 41~50)  
Improvement of abdominal tone, strengthening of the abdominal muscle development of firmer abdomen.

## 3. Contraindications and adverse reactions

### 3.1 Contraindication

- a) The long-term effects of chronic electrical stimulation to the human body are not known.
- b) Patients with cardiac pacemaker are disabled!
- c) Patients with suspected or confirmed have epilepsy, heart disease, pregnancy or menstruating women are prohibited.
- d) Treatment parts have the haemorrhage after acute trauma or fractures, disabling wound healing stage after the operation.
- e) Treatment is not sensitive to heat or electrical stimulation of the skin patients disabled.
- f) Disturbance of consciousness or disabled children.

- g) Treatment of swelling, skin infection or with skin disease, phlebitis patients disabled.
- h) With wet skin after shower, sweating and disabled when you go to sleep.
- i) Driving traffic tools or disable during movement.
- j) Metal allergies are disabled.

### 3.2 Adverse Reaction

- 1) You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin;
  - 2) You may experience headache and other painful sensations during or following the application of electrical stimulation near your eyes and to your head and face;
- You should stop using the device and should consult with your physician if you experience adverse reactions from the device.

## 4. Primary Structure of Product

The Transcutaneous Electrical Nerve Stimulators mainly consists of main unit, electrode patches and electrode wires. The primary structure of main unit as shown below:


Fig 1. Primary structure of main unit

## 5. Product Main Features

- The design is on the basis of ergonomics and is beautiful science, compact easy to carry;
- The LF electric pulse come into being an efficient composite energy field;
- Multi-function and multi-site remission,
- There are alternative functional modules, adjustable treatment time, and adjustable intensity;
- There is LCD display.

## 6. Product Requirements and Major Parametric Description

### 6.1 Product Power Requirement

- a) Power source: 3.7V/250mAh lithium battery
- b) Device safety classes: class II type BF

### 6.2 Product Main Technical Parametric description

- a) Impulse frequency: TENS:(20-100)Hz, EMS:(1-15)Hz, FITNESS:(2-16)Hz.
- b) Pulse width: TENS: 120 $\mu$ s, EMS: 200 $\mu$ s, FITNESS: 200 $\mu$ s.
- c) Impulse waveform: square wave
- d) The individual pulsing electric quantity when in the maximum output amplitude: >7 $\mu$ C.
- e) The maximum output energy of individual pulse:  $\leq$ 300mJ.
- f) Influence of output end open circuit and short circuit: it is able to support the influence of output end open and short circuit, and its performance will not be allowed to weaken.
- g) Adjustment of output amplitude: 0-16 levels.
- h) The timer of therapeutic equipment: 15 minutes, 30 minutes, 45 minutes,  $\pm$ 10% franchise.
- i) Safety classification: internal electric source class and type BF equipment
- j) Boundary dimension:  
Main unit: 122mm \* 55mm \* 15.7mm;  
Electrode Pads: 3 kinds  
EPAD-D01: 70 $\times$ 52mm EPAD-D02: 50 $\times$ 50mm EPAD-D03: 50 $\times$ 50mm
- k) Service life: the shelf life of subject device is 2 years; the use life of the small Electrode Pad is 80 times and shelf life is 2 years

### 6.3 Product Environmental Requirement

- a) Normal work environmental requirements:  
Environment temperature: +5 $^{\circ}$ C--+40 $^{\circ}$ C;  
Environment humidity: 15%-93%RH;
- Atmospheric environment conditions: 700hPa-1060hPa.
- b) Storage environment requirements:  
Environment temperature: -25 $^{\circ}$ C--+70 $^{\circ}$ C;  
Environment humidity: 0-93%RH; Atmospheric
- Environment conditions: 700hPa-1060hPa.
- c) Transport environment requirements:  
Environment temperature: -25 $^{\circ}$ C--+70 $^{\circ}$ C;  
Environment humidity: 0-93%RH; Atmospheric
- Environment conditions: 700hPa-1060hPa.

Note: It needs at least 30 minutes for device to return to normal operating temperature from the minimum/maximum storage temperature.

## 7. Method of application

### 7.1 Power supply operation

- 1) Charge the lithium battery by connecting the device to an adapter.
- 2) Connect the components. According to Fig 1, connect the main unit and electrode pads with the electrode wires.

#### 3) Turn on/off

Press the "on/off" key can turn on the main unit, the power light is lightened, the working routine is started and the LCD backlight of remote control is lightened.

Press the "On/Off" button again, the machine shut down and the screen is off.

### 7.2 Functional operation

- 1) Attach the electrode pads to the stimulate parts, if there is a perspiration in your body, it must be clean to use.
- 2) Press the "on/off" key can turn on the main unit. The device would be default mode of TENS, there are 20 fixed programs.
- 3) The use of the key " ^ "、" V " of program selection:  
" ^ "key as an option to Up select; " V " key as an option to Down select
- 4) The use of the key "+", "-" of intensity adjustment for channel 1 and channel 2.  
Each channel has 0-16 levels.  
"+" can enhance intensity; "-" can decrease intensity.
- 5) Mode switch  
Short press "F" key to switch three major modes of TENS, EMS, and FITNESS;
- 6) The adjustment of timing  
Long press "F" key for 1.5 seconds to set up the stimulation time of 15 minute, 30 minute, and 45 minute;
- 7) There are 20 fixed programs for TENS, 20 fixed programs for EMS and 10 fixed programs for FITNESS.
- 8) The output frequency of the product is displayed on the LCD.

### 7.3 The use of battery

- 1) The device powered by the lithium battery, which can rechargeable through the adapter.
- 2) When pressing the key, the brightness of LCD is distinct weakening or the stimulation intensity changes is very slightly when press the key of intensity adjustment, showing the electric quantity of battery is not enough and please charging in time.
- 3) If it is not used for one month or long term, please fully charge it first.
- 4) Do not use the device under the environment which more than 45 °C, otherwise it will affect the performance and life of the battery.
- 5) The discarding method of the post-batteries should be deal with according to the urban environmental protection.

## 8. Safety Precautions

- The manual said that the purpose of warning sign and legend is safe and proper using the product by yours, and prevents the harm to you and others.
- The warning sign and legend as well as their meaning is as follows:

Warning Sigh	Meaning
 using contraindication	It shows that it will appear the dangerous of casualties or serious injury in the error use.
 Warning	It shows that it will appear the possibility of casualties or serious injury in the error use.
 attention	It shows that it will appear the possibility of personal injury or damaged goods in the error use.

 attention	
<ul style="list-style-type: none"> <li>• The therapeutic apparatus cannot be used with the HF apparatus to avoid burns or damaged apparatus</li> <li>• If the patient himself uses the therapeutic apparatus and HF apparatus at the same time, the part of massage plate may cause burns on apparatus, it may also damage the apparatus;</li> <li>• If use the apparatus near (1 meter) the short wave or microwave Therapeutic apparatus, the output of apparatus may be unstable.</li> <li>• It will increase a danger of heart fibrillation by using electrode pads close to the chest.</li> <li>• Do not modify this equipment without authorization of the manufacturer</li> <li>• When need replacing the lithium batteries, please contact aftersales staff designated and authorized by the manufacturer directly, its accessories shall not be replaceable.</li> <li>• The Lithium battery is not replaceable by the user. Do not disassemble this device.</li> <li>• Shall not make personnel can no longer take care of themselves, infant, or not sensitive person use.</li> <li>• Simultaneous connection of a patient to a high frequency surgical ME equipment may result in burns at the site of the stimulator electrode pads and possible damage to the stimulator.</li> <li>• Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy ME equipment may produce instability in the stimulator output.</li> <li>• Application of electrode pads near the thorax may increase the risk of cardiac fibrillation.</li> <li>• Advice that stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrode pads placed on the chest and the upper back or crossing over the heart.</li> <li>• Prevent inhalation or accidental swallowing of small parts</li> </ul>	

<ul style="list-style-type: none"> <li>• Prevent sharp parts from damaging the product.</li> <li>• Do not use accessories not specified by the manufacturer.</li> </ul>	
 using contraindication	
<p>Pregnant women and women's menstrual period, person with sensitive skin, heart disease, abnormal blood pressure, malignant tumors, cerebrovascular patients, patients with acute disease or other person treated by doctor must consult a doctor before they can use this product.</p>	
<ul style="list-style-type: none"> <li>• (1) It is contraindicated for use by those who has a skin perceptual disturbance or is not sensitive about the heat.</li> <li>• (2) It is prohibited to use when bathing, sweating and sleeping.</li> <li>• (3) The patient with cerebral hemorrhage: It should be disabled in unsteady phase; the person who has sequela must be used under supervision of doctors.</li> <li>• (4) It is contraindicated for use by those who has purulent inflammation, acute blood poisoning and continuous hyperpyrexia.</li> <li>• (5) It is contraindicated for use by those who has acute cardiovascular and cerebrovascular diseases.</li> </ul> <p>Please stop using it immediately to the doctor to consult when unwell felt or skin complaint in the process of use.</p>	
 precautions	
<p>Please do not use in the parts which is nearest heart, head, eyes, front neck (especially the carotid artery), lower back, oral cavity or pudendum, skin disease.</p>	
<ul style="list-style-type: none"> <li>• (1) It should be shut down then using again when the apparatus is moved or changed the therapeutic parts in the using process, otherwise there will be a strong stimulation.</li> <li>• (2) It is not permitted to give children or people who were no ability to express consciousness to use.</li> <li>• (3) If you feel unwell because of using the product, please stop using it immediately to the doctor consult.</li> <li>• (4) Please pull up the power plugs when you are done or do not use the product.</li> <li>• (5) Do not used in conjunction with other medical electronic devices, such as cardiac pacemaker, artificial heart and lung and other</li> </ul>	

medical electronic devices, such as cardiac pacemaker, artificial heart and lung and other medical electronic devices with fuels life, electrocardiograph and other medical electronic devices, otherwise it will lead to danger.

(6) Do not use the product in the place where is direct sunlight, high heat,

- inflammable, electromagnetic radiation and humid.

(7) Do not disassemble, repair and transform the therapeutic apparatus, otherwise

- cause failure or get an electric shock accident.

(8) The therapeutic apparatus should be placed in the position which is easy to

- move the plugs to use, it is easy to move the plugs in emergency circumstances.

(9) Check the equipment before each use to avoid the exposed wires caused by accidental damages or other reasons.

- (10) Dust may affect the performance of the unit, please use a dry soft cloth to clean the device as needed. Please check whether the electrode is loose before each use, otherwise it may have adverse effects on performance or cause other problems.

## 9. Product maintenance

1) Before turning on the main unit, please check whether the battery is charged;

2) If the apparatus runs normally, show key invalidly after press. To turn it off after a few seconds and then opening, if still not normal, check whether it is broken for the apparatus;

3) If all operation of the apparatus are normal but without output (the user without feelings) after turn it on. Please check whether the electrode pads is directly contacted with skin, the contacting parts whether has hair, clothes and so on;

4) While cleaning, shutdown first, let the electrode wires stay away from the electric supply socket. You can wipe by using soft cloth or towel with a little water and after dry.

5) After each use, please clean electrode pads, you can wipe by using soft cloth or towel with a little water and after dry, when conducting electrode pads are particularly dirt, you can wipe by using soft cloth with a little medicinal alcohol (alcohol concentration as 75%).

6) The apparatus should be placed in the place which is dry, ventilated and good insulation use.

7) When the apparatus is moved, must be handled with care, avoiding the shake.

8) The apparatus should be checked the batteries condition, to check whether the output is good, if has abnormal conditions, you should pay attention to it and service it in time.

9) The apparatus' shell board should be pay attention to protection, avoiding the abrasion.

10) The non-professional cannot disassemble the host to avoid getting an electric shock or damaging the apparatus, an accident can happen because of it, all at your peril.

## 10. Product Scrap Processing

The therapeutic equipment belongs to medical device, if the complete machine is already aging and is not able to use, the process mode according to the local legislation carry through scrap processing.

## 12. Specifications

Model No.: UM-004

Product Name:

3 in 1 TENS EMS FITNESS Electrostimulator

Características eléctricas: 5V 180mA

## 11. Box Contents

1 x Main unit

4 x Electrode pads

1 x Charging Cable

2 x Electrode wires




1 x User Manual

NOTE: There are 3 kinds of replacement electrode pads in the below table. You can choose the most suitable, according to the electrode wires of your model

Model	Drawing	Size	Effective area of single piece (cm <sup>2</sup> )	Material
EPAD-D01		70x52 mm	33 cm <sup>2</sup>	Hydrogel
EPAD-D02		50x50 mm	25 cm <sup>2</sup>	Hydrogel
EPAD-D03		50x50 mm	25 cm <sup>2</sup>	Hydrogel

## 13. Meaning of the Symbols

Graphic symbol	Meaning
	Product batch number tag
	Product serial number tag
	Manufacturer
	Production date
	Precautions
	EU representative
	Warning against danger
	Applied part of type BF
	Class II equipment
	This product only for indoor use
	It means prohibition (the thing was not permitted to do)
	It means enforcement (the thing must be observed)
	Keep out of the sun
	Keep away from rain
	CE label and notified body code
	Contact its local authorities to determine the proper method of disposal of potentially bio hazardous parts and accessories.
IP22	The dustproofing and waterproofing grade, that is, to prevent the solid object is greater than 12mm to invade, the inclination of 15 degrees can prevent dripping invasion still not cause harmful effects.
	This symbol is used to direct the user to refer to documentation for additional information regarding the system use or description.

## 14.Executive Standards

The product conforms to the following standards and laws:

- IEC 60601-1:2005+A1: 2012 Medical electrical equipment-Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-11: 2015 Medical Electrical Equipment - Part 1-11: General Requirements For Basic Safety And Essential Performance - Collateral Standard: Requirements For Medical Electrical Equipment And Medical Electrical Systems Used In The Home Healthcare Environment
- IEC 60601-2-10: 2013 Medical electrical equipment - Part 2-10: Particular requirements for the safety of nerve and muscle stimulators
- IEC 60601-1-2:2014 Medical Electrical Equipment - Part 1-2: General Requirements For Basic Safety And Essential Performance - Collateral Standard: Electromagnetic Disturbances - Requirements And Tests

## Appendix EMC Declaration

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the EQUIPMENT should assure that it is used in such an environment.

The Transcutaneous Electrical Nerve Stimulators is suitable for use in a professional health care environment, not including areas where there are sensitive equipment or sources of intense electromagnetic disturbances, such as the RF shielded room of an imaging system magnetic resonance imaging, in operating rooms near active AF surgical equipment, electrophysiology laboratories, armored rooms or areas where short wave therapy equipment is used.

- Do not use the system around strong electric field, electromagnetic field (e.g. MRI scan room) and mobile wireless communication devices. Using the device in an improper environment may cause malfunction or damage.
- The compliance with EMC and EMI regulation cannot be guaranteed by the use of modified cables or those which does not comply with the same standards under what the equipment was validated.
- The system must not be used adjacent or supported by other equipment. The recommendations of this manual must be followed.
- Do not use accessories, transducers, internal parts of components and other cables other than those previously specified by the manufacturer. This may result in increased emission or decreased electromagnetic immunity and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm to any part of the ultrasound system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- To maintain basic safety in relation to electromagnetic disturbances during the expected service life, always use the system in the specified electromagnetic environment and follow the maintenance recommendation described in this manual.

Danger

The following tables provide information on compliance of the equipment according to the standard EN 60601-1- 2:2015.

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Emissions Test	Compliance	Electromagnetic Environment and Guidance
RF emissions CISPR 11	Grupo 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2- Compliance standards

Phenomenon	Basic Standard of	Immunity Test	Level of
Electrostatic discharge	IEC 61000-4-2	±8 KV contact ±2 KV, ±4 KV, ±8 KV,±15KV air	±8 KV contact ± 2 KV, ± 4 KV, ± 8 KV, ±15KV air
Radiated RF EM fields <sup>1</sup>	IEC 61000-4-3	3V/m 80 MHz-2.7 GHz 80% AM at 1 KHz	3V/m 80 MHz-2.7 GHz 80% AM at 1 KHz
Proximity fields from RF wireless communication equipment	IEC 61000-4-3	See table	See table
Electrical Fast/ Transients bursts	IEC 61000-4-4	±1 KV 100 KHz repetition frequency	±1 KV 100 KHz repetition frequency
Conducted disturbances induced by RF fields.	IEC 61000-4-6	3V 0.15 MHz-80 MHz 6 Vm in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1KHz	3V 0.15 MHz-80 MHz 6 Vm in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1KHz
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz

**Table 3- Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipments**

Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM $\pm 5$ KHz deviation 1KHz sine	2	0.3	28
710	704-787	LTE 13, 17 Band	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE 5 Band	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700- 1990	GSM 1800 , CDMA 1900, GSM 1900, DECT,LTE 1, 3, 4 , 25 Band, UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400- 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE 7 Band	Pulse modulation 217 Hz	2	0.3	28
5240	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

### OZU Warranty

OZU guarantees this product for a period of 24 MONTHS from the date of purchase. Defects in materials and workmanship will be repaired or replaced by OZU, provided valid proof of purchase is presented within the specified period. The use of unauthorized replacement parts will void this warranty. The user must strictly comply with all the instructions in the instruction manual and must refrain from any action or use that is considered inappropriate or that is warned in it. If the product is used in a country where it should not be used, voltage conversion may be required. OZU is not responsible for power failures due to incorrect voltage supply to the device.

### Warranty conditions

The warranty begins on the date of purchase and expires at the end of the period indicated above. A valid proof of purchase is required. OZU's warranty applies if the product has been handled properly for its intended use and used in accordance with the operating instructions at all times.

The warranty does not cover losses caused by nature, including but not limited to data loss or business loss. Some warranty exclusions apply.

### Warranty Terms

The warranty begins on the day of purchase and expires at the end of the period indicated above. In the event that the merchant from which it was purchased is not an authorized OZU channel (for example, on online auction sites), the warranty period is considered to start three months after the date of manufacture indicated on the product or on the batch or serial number. Should a defect occur due to faulty workmanship or materials within the warranty period, OZU will repair it free of charge. Where repair is impossible or deemed uneconomical, OZU may agree to replace the product, at OZU's discretion. The OZU warranty will apply as long as the product has been used correctly according to the user manual provided and the purpose it was produced for. All expenses derived from transportation to and from the service center will be borne by the customer.

### Exclusions

The warranty does not cover losses due to natural causes, including but not limited to loss of data or business.

The causes of warranty cancellation are as follows:

- The purchase document was modified or altered in any way.
- The serial number, batch, or model of the product is modified, altered, removed, or made illegible.
- A repair, modification, or alteration has been carried out by unauthorized personnel or unauthorized technical service.
- The product was used for commercial purposes or other than domestic use.
- The defect was caused by abuse or misuse of the product or by environmental conditions that do not agree with what is recommended for its operation.
- The defect was caused by connection to peripheral equipment, additional equipment, or accessories other than those recommended by OZU.
- The unit has been damaged, including, but not limited to damage caused by pets, lightning, abnormal voltage, water or fire, natural disasters, or transportation

accidents.

- Non-functional parts are considered consumable parts by their nature and are not included in the warranty.
- The product does not work correctly because it was not designed, manufactured, approved, and/or authorized for use in the country where it was used.
- Excesses or drops in electrical voltage that imply use in abnormal conditions, as well as defects in the electrical installation.
- Adulteration, absence, or amendment of the serial number or date of manufacture or of the label that contains the model, as appropriate.
- Damage caused by impacts, falls, improper transport and storage, or improper treatment, as well as damage due to tremors, earthquakes, tidal waves, floods, lightning, humidity, excessive dust environments, etc.
- The operation of the product in conditions not prescribed in the user manual.
- Technical defects resulting from problems in the incorrect installation of the product.
- If applicable, the use of products, equipment, systems, spare parts, accessories, applications, external connections, installations, or connectors not supplied or recommended in the user manual.
- OZU brand products that are not imported by Ventis SpA for the Chilean market, by Ventis Perú SAC for the Peruvian market, by Ventis Colombia SAS for the Colombian market, or by Ventis Retail S.A.P.I. de C.V. for the Mexican market.

## IMPORTANT

BEFORE USING THE PRODUCT, READ THE INSTRUCTIONS AND KEEP THEM FOR FUTURE REFERENCE.

## SAFETY MEASURES

When using electrical products, especially when children are present, basic safety precautions in the manual should be followed.

## Customer service

We recommend that you carefully read the product's operating instructions before contacting OZU or the vendor. To obtain service within the warranty period, please contact the corresponding Authorized Technical Service.

In order to assist you more efficiently, please have your purchase receipt or invoice ready when contacting OZU, indicating the date of purchase, vendor, model, and serial number, production batch, or production date as indicated on the product.

## Chile

+56 22 840 5567

WhatsApp: +56 9 92618728

## Peru

0800 56185

Mexico

+52 55 6826 9609

WhatsApp: +52 55 1591 9179

Colombia

+57 320 999 6287

Telephone service hours: Monday to Friday from 9:00 a.m. to 6:00 p.m.

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CHILE: Ventis SpA. Camino Lo Boza 107 A09, Pudahuel, Santiago, Chile. RUT 76229287-4
COLOMBIA: Ventis Colombia SAS. CR 10 N°96 25 OF 612. NIT 901375385